

HandyLab 7series THE NEW MOBILE MEMOSENS® MEASURING DEVICES



a **xylem** brand

20200 STATE BURNESSES ۵ 50 00 (III)

1-1-1

Copyright Note: MEMOSENS® and MemoSuite® are registered Wordmarks of Endress+Hauser Conducta GmbH + Co. KG, Gerlingen. The Logo MEMO SENS is a registered Trademark of Endress+Hauser Conducta GmbH + Co. KG, Gerlingen and Knick Elektronische Messgeräte GmbH & Co. KG, Berlin

Content

| The new HandyLabs | Page | 4 |
|--|------|----|
| Model Selection Table | Page | 7 |
| HandyLab 700 | Page | 8 |
| HandyLab 700 - Specifications | Page | 10 |
| HandyLab 750 | Page | 12 |
| HandyLab 750 - Specifications | Page | 14 |
| HandyLab 750EX | Page | 16 |
| HandyLab 750EX - Specifications | Page | 18 |
| HandyLab 780 | Page | 20 |
| HandyLab 780 - Specifications | Page | 22 |
| PC-Software - HandyLab Pilot | Page | 24 |
| Order Information for HandyLab and Accessories | Page | 26 |
| Memosens® Electrodes | page | 28 |
| Order Information Memosens®-Electrodes | page | 30 |

The new mobile pH measuring devices by SI Analytics with Memosens[®] technology offers increased safety and a userfriendly interface

MEMOSENS[®] - superior engineering

When measuring pH, both the sensor and electrical connection to the measuring device are critical for providing quality process data. TTraditional analog sensors contain high interior resistors that require a high-ohm connection to the measuring device. The presence of moisture on the sensor contacts can cause inaccurate readings and even sensor failure. This risk is eliminated by Memosens technology. The measured value is processed in the sensor head, digitalized, and then transmitted to the measuring device without metal contacts; there is no interference.

SI ADAINTICS

by SI Analytics

MEMOSENS[®] Electrodes

HandyLab 780 with Memosens®-Elektrode

4



Sensor and PC connections as well as electrode storage.

| Clock | |
|--|--|
| on/off | |
| | |
| Managara seate i 2 2 2 j information Zarrick <u>ensity</u> - <u>()</u> - | |

This is how you secure your data

With conventional analog systems the sensor must be calibrated with every replacement. This is due to the calibration data being stored only in the measuring device. The Memosens concept solves this problem by saving the calibration data in the sensor head. When connected to the measuring device, they automatically identify themselves with their serial number and part number and transmit their calibration data to the measuring device. There is no additional calibration required to ensure an accurate measurement.

- Increased safety due to digital signal transmission. External interferences are eliminated, calibration data is transferred securely, and sensor data is transmitted easily.
- The sensor data is stored in the sensor.
- Predictive maintenance thanks to the possibility to track the sensor's past performance.
- Submersible thanks to its hermetically sealed plug-in head.
- Memosens[®] is an open system, which means that it is supported by several manufacturers and is a defacto standard.

HandyLab 7series

The new mobile Measuring Devices by SI Analytics



Selection Table HandyLab 7series

| Application | HandyLab 700 | HandyLab 750 | HandyLab 750EX | HandyLab 780 |
|------------------------------------|--------------|--------------|----------------|---------------------------|
| Memosens® pH and ORP (Redox) | | | | |
| Memosens [®] Conductivity | | | | |
| Memosens® Oxygen | | | | |
| Analog pH and ORP (Redox) | | | | |
| Temperature | | | | |
| Ex-Zone 0/1 | _ | | | - |
| PC Program HandyLab Pilot | _ | | | |
| Micro USB-B | _ | | | |
| Data logger (values) | - | 5,000 | 5,000 | 10,000 |
| Li-Ion battery | - | | _ | |
| Display | LCD segment | LCD segment | LCD segment | OVGA-TFT Color graphic |
| Multiple languages | _ | _ | _ | |
| Help functions | _ | _ | _ | |

HandyLab 7series - offers increased safety and a userfriendly interface:

- Brings the advantages of digital Memosens[®] technology, offering increased safety and a userfriendly interface to laboratories, technical colleges, field and process.
- Durable and chemical resistant housing
- Has passed the standardized drop test from 1 m height on cement
- The protective cap protects the device.
- The HandyLabs can be hung up on a hook integrated into the housing.
- By flipping the protective flap and holding hooks, these HandyLab devices can either be suspended or used as bench top units.
- The integrated electrode vessel protects the sensors from drying out or being damaged and is detachable from the housing.
- Simultaneous pH/mV and temperature display
- Can be connected to Memosens[®] as well as analog sensors
- Protection class IP 67/ IP 66 (splash water protected)





HandyLab 700

Entering the world of mobile Memosens®





Durable design for daily measurements

- Contrast-rich and scratch-resistant clear glass LCD display
- Self-explanatory operation with one clear text line
- Display of the sensor status at one glance with Sensoface
- Calibration of up to 3 points with automatic buffer detection (from 10 buffers)
- Manual calibration by indicating random buffer values
- Real-time clock and display of battery charging status
- Extended operating time of more than 1,000 hours with standard batteries (4 x AA)







HandyLab 700 - Specifications

| Measuring parameters | | | |
|--|--|---|--|
| Memosens pH (also called ISFET) | Connection: Female connecto Display areas ¹⁾ | or M8 (4-way) for Memosen: -2.000 +16.000 pH -2000 +2000 mV -50 +250 °C | s laboratory cable |
| Memosens Redox | Connection: Female connecto Display areas ¹⁾ Sensor adaptation ^{*)} Perm. calibration range | or M8 (4-way) for Memosen: -2000 +2000 mV -50 +250 °C ORP (Redox) calibration (ze ΔmV (Offset) -700 +700 | s laboratory cable ero point offset) mV |
| pH/mV (analog) | Connection Measuring Range pH Decimal places [*]) Input Resistance Input current Measuring cycle Operation measurement deviation ^{2, 3, 4}) Measuring range mV Measuring cycle Operation measurement deviation ^{2, 3, 4}) | pH female connector DIN 19 262 (13/4 mm) -2 16 2 or 3 1 x $10^{12} \Omega$ (0 35 °C) 1 x $10^{-12} A$ (with RT, double approx. 1s < 0.01 pH, TK < 0.001 pH/I -1300 +1300 approx. 1s < 0.1 % v. M. + 0.3 mV TK < | e up every 10 K) < < 0.03 mV/K |
| Temperature | Connection: 2 x Ø 4 mm for in Measuring ranges Measuring cycle Operation measurement deviation ^{2, 3, 4)} | ntegrated or separate temp NTC 30 kΩ Pt 1000: approx. 1s < 0.2 K (Tamb = 23 °C); TK | erature sensor -20 +120 °C -40 +250 °C < 25 ppm/K |
| pH Calibration | | | |
| Operating Modes ^{*)} AutoCal Buffer Sets ^{*)} | AutoCal - Calibration with aut Manual - manual calibration v Data entry - data entry of zero Knick CaliMat | comatic buffer detection while entering individual buf point and incline Ciba (94) | fer values Hamilton |
| | NIST Technical | НАСН | Mettler-Toledo |
| | NIST-Standard DIN 19267 | SI-Analytics techn. buffer Reagecon | |
| Perm. calibration range | Zero point | 68рН | |
| | For ISFET: | -750 +750 mV work poi | nt (asymmetrical) |
| | Slope | approx. 74 104 % | |
| Calibration Timer ^{*)} | Sample interval 1 99 days, | can be switched off | |
| Display - Operation | | | |
| Sensoface | Provides information regardir incline, setting time, calibratic | ng the sensor status, evaluat on interval status display (frie | ion of zero point/ endly, neutral, sad). |
| Display | LCD STN 7-segment display w | with three lines and symbols | ; |
| Status displays | For battery status | | |
| Notes | Sand timer | | |
| Keyboard | [on/off], [cal], [meas], [set], [|], [▼], [clock] | |

| Diagnostic Functions | Sensor data (Memosens only) Manufacturer, sensor type, serial number, operating time | |
|--------------------------------------|---|--|
| | Calibration data | Calibration date, zero point and slope |
| | Device self-test | Automatic memory test (FLASH, EEPROM, RAM) |
| | Machine data | Device type, software version, hardware version |
| Data retention | Parameters, calibration data > | • 10 years |
| Climate - rated operating conditions | Ambient temperature Transport/storage temperature Relative humidity | -10 +55 °C -25 +70 °C 0 95 %, short-term dew permitted |
| Energy supply | Auxiliary energy | Batteries 4x AA (Mignon) alkaline or lithium |
| | Operating time | approx. 1,000 h (Alkaline) |
| Casing | material | PA12 GF30 + TPE |
| | Type of protection | IP66/67 with pressure compensation |
| | Dimensions | approx. 132 x 156 x 30 mm |
| | Weight | approx. 500 g |
| Certificates - testing mark - d | levice safety | |
| EMV | DIN EN 61326-1 (General Rec | quirements) |
| | Interference Emission | Class B (residential) |
| | Interference resistance | Industry Branch |
| | DIN EN 61326-2-3 (Special Requirements for Pressure Transducers) | |
| RoHS Conformity | As per directive 2011/65/EU | |

- *) parameterizable
 1) Measuring ranges depend on Memosens sensor
 2) As per DIN EN 60746-1, with rated operating conditions
 3) ± 1 digit
 4) Plus sensor errors

11

HandyLab 750

Portable Memosens[®] world with data storage and USB interface



- Store up to 5,000 data points
- HandyLab Pilot software for evaluation of data from digital Memosens[®] sensors.
- HandyLab Pilot software (included with delivery), allows easy management and evaluation of stored data and can be used for device configuration.
- Contrast-rich and scratch-resistant clear glass display.
- Self-explanatory operation with one clear text line
- Display of the sensor status at one glance with Sensoface.
- Calibration of up to 3 points with automatic buffer detection.
- Manual calibration by indicating random buffer values
- Real-time clock and display of battery charging status.
- Extended operating time of more than 1,000 hours with standard batteries (4 x AA) or use of a Li-Ion accumulator, even with very high or low operating temperatures. The battery is charged via the USB interface.







HandyLab 750 - Specifications

| Measuring parameters | | | |
|--|--|--|---|
| Memosens pH (also called | Connection: Female connector M8 (4-way) for Memosens laboratory cable | | aboratory cable |
| | Display areas ¹⁾ | -2,000 +16,000 pH -2000 +2000 mV -50 +250 °C | |
| Memosens Redox | Connection: Female connector Display areas ¹⁾ Sensor adaptation ^{*)} Perm. calibration range | M8 (4-way) for Memosens® -2000 +2000 mV -50 +250 °C ORP (Redox) calibration (ze AmV (Offset) -700 +700 | lab cable ero point offset) mV |
| pH/mV (analog) | Connection: pH female connect Measuring Range pH Decimal places ^{*)} Input Resistance Input current Measuring cycle Operation measurement deviation ^{2, 3, 4)} Measuring range mV Measuring cycle Operation measurement deviation ^{2, 3, 4)} | tor DIN 19 262 (13/4 mm) -2 16 2 or 3 1 x $10^{12} \Omega$ (0 35 °C) 1 x $10^{-12 A}$ (for RT, doubles every 10 K approx. 1s < 0.01 pH, TK < 0.001 pH/f -1300 +1300 mV approx. 1s < 0.1 % v. M. + 0.3 mV; TK | .) K < 0.03 mV/K |
| Temperature | Connection: 2 x Ø 4 mm for inte Measuring ranges Measuring cycle Operation measurement deviation ^{2, 3, 4}) | egrated or separate temper NTC 30 kΩ Pt 1000: approx. 1s < 0.2 K (Tamb = 23 °C); TK | ature sensor -20 +120 °C -40 +250 °C < 25 ppm/K |
| pH Calibration Operating Modes ^{*)} AutoCal Buffer Sets*) | AutoCal Manually Data entry Knick CaliMat NIST Technical | Calibration with automatic Manual calibration while er buffer values Data entry of zero point an Ciba (94) HACH | buffer detection ntering individual d incline Hamilton Mettler-Toledo |
| Perm. calibration range | NIST-Standard DIN 19267 Zero point For ISFET: Slope | SI Analytics techn. buffer Reagecon 6 8 pH -750 +750 mV; work po approx. 74 104 % | int (asymmetric) |
| Calibration Timer" | Sample interval 1 99 days, ca | an be switched off | |
| Display - Operation Sensoface Display Status displays | provides information regarding Evaluation of zero point/slope, Status display (friendly, neutral, LCD STN 7-segment display wi For battery status, logger | the status of the sensor setting time, calibration inte sad) ith three lines and symbols | rval |
| Notes | Sand timer | | |
| Keyboard | [on/off], [cal], [meas], [set], [▲], | , [▼], [STO], [RCL], [clock] | |

| Diagnostic Functions | Sensor data (Memosens only) | Manufacturer, sensor type, serial number, operating time |
|--|--|--|
| | Calibration data | Calibration date, zero point and slope |
| | Device self-test | Automatic memory test (FLASH, EEPROM, RAM) |
| | Machine data | Device type, software version, hardware version |
| Data retention Data transmission Data logger | Parameters, calibration data > 1 1x Micro USB-B for data transmi 5.000 memory positions | 0 years ission to the PC |
| | Recording | controlled manually, by intervals or events |
| Calibration data logger MemoLog | up to 100 Memosens calibration | n logs can be saved |
| (Memosens only) | - Recording can be displayed o - directly readable via MemoSu zero point, incline, calibration d | n the display ite® (USB): Manufacturer, sensor type, serial no., ate |
| Communication | USB 2.0 | |
| | Profile Use | HID, installation without driver Data exchange and configuration via the software HandyLab Pilot |
| Climate - rated operating con | ditions | |
| | Ambient temperature | -10 +55 °C |
| | Transport/storage temperature | -25 +70 °C |
| | Relative humidity | 0 95 %, short-term dew permitted |
| Energy supply | | |
| | Auxiliary energy | Batteries 4x AA (Mignon), 4x Akku NiMH or 1x Li-Ion battery pack, chargeable via USB |
| | Operating time | approx. 1,000 h (Alkaline) |
| Casing | | |
| | material | PA12 GF30 + TPE |
| | Type of protection | IP66/67 with pressure compensation |
| | | approx. 132 x 156 x 30 mm |
| Cartificates testing mande de | Weight | approx. 500 g |
| EMV | DIN EN 61326-1 (General Requ | irements) |
| | Interference Emission | Class B (residential) |
| | Interference resistance DIN EN 61326-2-3 (Special Req | Industry Branch uirements for Pressure Transducers) |
| RoHS Conformity | As per directive 2011/65/EU | |

*) parameterizable
1) Measuring ranges depend on Memosens sensor
2) As per DIN EN 60746-1, with rated operating conditions
3) ± 1 digit
4) Plus sensor errors

HandyLab 750EX

Portable Memosens[®] world with data storage and USB interface for the use in potentially explosive areas



- For use in potentially explosive areas up to zone 0/1.
- Store up to 5,000 data points.
- Micro-USB connection to communicate with the software HandyLab Pilot software for evaluation of data from digital Memosens[®] sensors.
- HandyLab Pilot software (included with delivery), allows easy management and evaluation of stored data and can be used for device configuration.
- Contrast-rich and scratch-resistant clear glass display.
- Self-explanatory operation with one clear text line.
- Display of the sensor status at one glance with Sensoface.
- Calibration of up to 3 points with automatic buffer detection.
- Manual calibration by indicating random buffer values.
- Real-time clock and display of battery charging status.
- Extended operating time of more than 1,000 hours with standard batteries (4 x AA).







HandyLab 750EX - Specifications

| Measuring parameters | | | |
|---|--|---|--|
| Memosens pH (also called | Connection: Female connector M8 (4-way) for Memosens laboratory cable | | |
| | Display areas ¹⁾ | -2.000 +16.000 pH -2,000 +2,000 mV -50 +250 °C | |
| Memosens Redox | Connection: Female connecto Display areas ¹⁾ Sensor adaptation ^{*)} Perm. calibration range | or M8 (4-way) for Memosens -2,000 +2,000 mV -50 +250 °C ORP (Redox) calibration (ze AmV (Offset) -700 +700 | s laboratory cable ero point offset) mV |
| pH/mV (analog) | Connection: pH female conne Measuring Range pH Decimal places ^{*)} Input Resistance Input current Measuring cycle Operation measurement deviation ^{2, 3, 4)} Measuring range mV Measuring cycle Operation measurement deviation ^{2, 3, 4)} | ector DIN 19 262 (13/4 mm) -2 16 2 or 3 1 x 10 ¹² Ω (0 35 °C) 1 x 10 ⁻¹² A (with RT, double approx. 1s < 0.01 pH, TK < 0.001 pH/I -1,300 +1,300 mV approx. 1s < 0.1 % v. M. + 0.3 mV TK < 0.03 mV/K | e up every 10 K) K |
| Temperature | Connection: 2 x Ø 4 mm for in Measuring ranges Measuring cycle Operation measurement deviation ^{2,3,4)} | ntegrated or separate tempo NTC 30 kΩ Pt 1000: approx. 1s < 0.2 K (Tamb = 23 °C); TK | erature sensor -20 +120 °C -40 +250 °C < 25 ppm/K |
| pH Calibration Operating Modes ^{*)} | AutoCal Manually Data entry | Calibration with automatic Manual calibration while er values Data entry of zero point an | buffer detection ntering individual buffer d incline |
| AutoCal Buffer Sets ^{*)} | Knick CaliMat NIST Technical NIST-Standard DIN 19267 | Ciba (94) HACH SI Analytics techn. buffer Reagecon | Hamilton Mettler-Toledo |
| Perm. calibration range | Zero point For ISFET: Slope | 6 8 pH -750 +750 mV; work po approx. 74 104 % | int (asymmetric) |
| Calibration Timer ^{*)} | Sample interval 1 99 days, | can be switched off | |
| Display - Operation | | | |
| Sensoface | Provides information regardi setting time, calibration inter | ng the sensor status, evalua val status display (friendly, n | tion of zero point/slope, eutral, sad). |
| Display | LCD STN 7-segment display | with three lines and symbols | S |
| Status displays | for battery status, loager | | |
| Notes | Sand timer | | |
| Keyboard | [on/off], [cal], [meas], [set], [|], [▼], [STO], [RCL], [clock] | |

| Diagnostic Functions | Sensor data (Memosens only) | Manufacturer, sensor type, serial number, operating time |
|------------------------------------|---|--|
| | Calibration data | Calibration date, zero point and slope |
| | Device self-test | Automatic memory test (FLASH, EEPROM, RAM) |
| | Machine data | Device type, software version, hardware version |
| Data retention | Parameters, calibration data > | 10 years |
| Data transmission | 1x Micro USB-B for data transı | mission to the PC |
| Data logger | 5,000 memory positions | |
| | Recording | controlled manually, by intervals or events |
| Calibration data logger MemoLog | up to 100 Memosens calibrati | on logs can be saved |
| (Memosens only) | - Recording can be displayed | on the display |
| Communication | - directly readable via MemoS point, incline, calibration date | Suite (USB): Manufacturer, sensor type, serial no., zero |
| Communication | Profile | HID, installation without driver |
| | Proper | Data exchange and configuration via the software HandyLab Pilot |
| Climate - rated operating c | onditions | |
| | Ambient temperature | $-10 \degree C \le Ta \le +40 \degree C$ T4 $-10 \degree C \le Ta \le +50 \degree C$ T3 |
| | Transport/storage temperature | -25 +70 °C |
| | Relative humidity | 0 95 %, short-term dew permitted |
| Energy supply | | |
| | Auxiliary energy | Batteries 4x AA (Mignon) |
| | Operating time | approx. 1,000 h (Alkaline) |
| Casing | | |
| | material | PA12 GF30 + TPE |
| | lype of protection | IP66/6/ with pressure compensation |
| | Dimensions | approx. 132 x 156 x 30 mm |
| | Weight | approx. 500 g |
| Certificates - testing mark - | device safety | |
| ENIV | DIN EN 61320-1 (General Rec | (lass B (residential) |
| | Interference resistance | Industry Branch |
| | DIN FN 61326-2-3 (Special Re | and the provide the provide the provided and the provided |
| Explosion protection | Furope | ATEX II 1 G Ex ia IIC T4/T3 Ga |
| RoHS Conformity | As per directive 2011/65/EU | |

- *) parameterizable
 1) Measuring ranges depend on Memosens sensor
 2) As per DIN EN 60746-1, with rated operating conditions
 3) ± 1 digit
 4) Plus sensor errors

HandyLab 780

Portable multi-parameter Memosens® world



Durable design for everyday use measurements

- With Memosens[®] sensor measurement of pH, conductivity and dissolved oxygen by simple sensor exchange
- Automatic detection of the measuring parameter with Memosens[®] sensors.
- Automatic compensation of the ambient pressure for oxygen measurement
- Storage of up to 10,000 data points.
- Micro-USB connection to communicate with the HandyLab Pilot software for evaluation of data from digital Memosens[®] sensors.
- HandyLab Pilot software (included with delivery), allows for easy management and evaluation of stored data and can be used for device configuration.
- Contrast-rich and scratch-resistant clear glass display.
- Self-explanatory operation with extensive information and help text.
- Display of the sensor status at one glance with Sensoface.



- Calibration of up to 3 points with automatic buffer detection.
- Manual calibration by indicating random buffer values.
- Real-time clock and display of battery charging status.
- Extended operating time of more than 1,000 hours with standard batteries (4 x AA) or Li-Ion accumulators, even with very high or low operating temperatures. The batteries are charged via the USB interface.

Menüauswahi

Inton

Zurück

~ 0

рн 6.76

Menü

SI Analytics

25

Wert h

HandyLab 780 - Specifications

| Measuring parameters | | | |
|--|--|--|---|
| Memosens pH (also called ISFET) | Connection: Female connector | M8 (4-way) for Memosens la | boratory cable |
| | Display areas ¹⁾ | -2.000 +16.000 pH -2000 +2000 mV | |
| | | Temperature -50 +250 °C | |
| Memosens Redox | Connection: Female connector Display areas ¹⁾ Temperature | M8 (4-way) for Memosens la mV -2000 +2000 mV -50 +250 °C | boratory cable |
| | Sensor adaptation ^{*)} Perm. calibration range | ORP (Redox) calibration (zer ΔmV (Offset) -700 +700 r | ro point offset) nV |
| pH/mV (analog) | Connection: pH female connec | tor DIN 19 262 (13/4 mm) | |
| | Measuring Range pH | -216 | |
| | Decimal places | 2 or 3 1 $\times 10^{12} \text{ O} (0 - 25 \text{ °C})$ | |
| | Input current | $1 \times 10^{-12} \Lambda$ (with RT double | $\mu \rho $ $\rho \nu \rho \eta (10 \text{ K})$ |
| | Measuring cycle | approx 1s | |
| | Operation measurement deviation ^{2, 3, 4)} | < 0.01 pH, TK < 0.001 pH/K | |
| | Measuring range mV | -1300 +1300 | |
| | Measuring cycle | approx. 1s | |
| | Operation measurement deviation ^{2, 3, 4)} | < 0.1 % v. M. + 0.3 mV TK < | 0.03 mV/K |
| Temperature | Connection: 2 x Ø 4 mm for inte | egrated or separate tempera | ture sensor |
| | Measuring ranges | NTC 30 kΩ -20 +120 °C | |
| | | Pt 1000 -40 +250 °C | |
| | Measuring cycle | approx. 1s | 0.5 |
| | Operation measurement deviation ^{2, 3, 4)} | < 0.2 K (lamb = 23 °C); IK < | < 25 ppm/K |
| Memosens Conductivity Connection: Female connector M8, 4-way for Memos | | M8, 4-way for Memosens® la | aboratory cable alternative |
| - | Measuring cycle | approx. İs | - |
| | Temperature compensation: | linear 0 20 %/K, reference nLF: 0 120 °C | e temperature adjustable |
| | | NaCl HCl (pure water with traces) | |
| | | NH_3 (pure water with traces NaOH (pure water with trace) |) es) |
| | Display resolution ⁵⁾ (autorangin | iq) | |
| | conductivity | 0.01 μS/cm (c < 0.5 cm-1) | |
| | - | $0.01 \mu\text{S/cm} (c = 0.05 \dots 0.2)$ | cm-1) |
| | | 0.1 µS/cm (c > 0.2 cm-1) | |
| | spec. resistance | 00.00 99.99 MΩ ∙ cm | |
| | Salinity | 0.0 45.0 g/kg (0 30 °C) |) |
| | TDS | 0 1999 mg/l (10 40 °C |) |
| | Concentration | 0.00 9.99 wt. % | |
| | Concentration definition | | 2) |
| | | $0.00 \dots 9.99$ wt. % $(0 \dots 60^{-6})$ | ~) °C) |
| | | $0.00 \dots 9.99$ wt % (0 100 ° | |
| | H-SO. | 0.00 9.99 wt % (-17 11) | |
| | HNO ₂ | 0.00 9.99 wt. % (-17 50 | °C) |
| | Sensor adaptation: | | _, |
| | Cell constant | Input of the cell constant wi | th simultaneous display of |
| | | the conductivity value and t | he temperature |
| | Enter solution | Input of the conductivity of the simultaneous display of the | the calibration solution with cell constant and the |
| | Auto | Automatic dotormination of | the cell constant with KCI |
| | Auto | solution or NaCl solution | |
| Memosens Oxygen | Connection: Female connector | M8, 4-way for Memosens lab | poratory cable |
| | Display areas ¹⁾ | Saturation | 0.000 200.0 % |
| | | Concentration | 000 μg/I 20.00 mg/l |
| | Manauring range to manareture 1 | Partial pressure | 0.0 1000 mbar |
| | Sensor adaptation: Automatic | alibration air and moisture a | diustable zoro point |
| | calibration | | ajustable, zero politi |

| pH Calibration | | | |
|---------------------------------|---|--|------------------------------------|
| Operating Modes ^{*)} | AutoCal | Calibration with automatic b | uffer detection |
| | Manually | Manual calibration while ent | ering individual buffer |
| | _ | values | |
| + | Data entry | Data entry of zero point and | slope |
| AutoCal Buffer Sets" | Knick CaliMat | Ciba (94) | Hamilton |
| | NIST Technical | HACH | Mettler-Toledo |
| | NIST-Standard | SI Analytics techn. buffer | |
| | DIN 19267 | Reagecon | |
| Perm. calibration range | Zero point | 68 H | |
| | For ISEET | -750 + 750 mV work point | (asymmetrical) |
| | Slope | approx 74 104 % | |
| Calibration Timer ^{*)} | Sample interval 1 99 days, ca | an be switched off | |
| Display - Operation | , | | |
| Device operation | clear menu guide with graphic instructions | symbols and detailed easy-to | o-understand operating |
| Languages | German, English, French, Spani | sh, Italian, Portuguese | |
| Sensoface | Provides information regarding setting time, calibration interval | the sensor status, evaluation status display (friendly, neutr | of zero point/slope, ral, sad). |
| Display | QVGA TFT display with white ba | acklight | |
| Status displays | for battery status, logger | | |
| Notes | Sand timer | | |
| Keyboard | [on/off], [meas], [enter], [▲], [▼ |], [◀], [▶], 2 Softkeys with cor | ntext-related population |
| Diagnostic Functions | Sensor data (Memosens only): | Manufacturer, sensor type, so operating time | erial number, wear, |
| | Calibration data: | Calibration date, zero point a | and incline |
| | Device self-test: | Automatic memory test (FLA | SH, EEPROM, RAM) |
| | Machine data: | Device type, software versio | n, hardware version |
| Data retention | Parameters, calibration data > 1 | 0 years | |
| Data transmission | 1x Micro USB-B for data transm | ission to the PC | |
| Data logger | 10,000 memory positions | | |
| | Recording: | controlled manually, by inter | rvals or events |
| Calibration data logger | up to 100 Memosens calibratio | n logs can be saved | |
| MemoLog | | | |
| (Memosens only) | - Recording can be displayed o | n the display | |
| | - directly readable via wemosu | Ite (USB): | ration data |
| Communication | lice 2.0 | ai no., zero point, siope, calibr | ation date |
| Communication | USD 2.0 Profile | HID installation without driv | or |
| | Proper | Data exchange and configur | ei ation via the software |
| | Порег | Handyl ab Pilot | |
| Climate - rated operating co | nditions | hanay2ab i not | |
| 1 5 | Ambient temperature | -10 +55 °C | |
| | Transport/storage temperature | -25 +70 °C | |
| | Relative humidity | 0 95 %, short-term dew pe | ermitted |
| Energy supply | - | | |
| | Auxiliary energy | Batteries 4x AA (Mignon), 4x | Akku NiMH or 1x Li-Ion |
| | | battery pack, chargeable via | USB |
| | Operating time | approx. 500 h (Alkaline) | |
| Casing | | | |
| | material | PA12 GF30 + TPE | |
| | lype of protection | IP66/6/ with pressure comp | ensation |
| | Dimensions | approx. 132 x 156 x 30 mm | |
| | | approx. 500 g | |
| Certificates - testing mark - o | | ize ne e nte) | |
| | UIN EN 01320-1 (General Requ | Clear D (marial anti-1) | |
| | Interference Emission | Class B (residential) | |
| | DIN EN 61226 2 2 (Special Devial | mousily dranch | ducorc) |
| Polls Conformity | As par directive 2011//E/EU | unements for Fressure fransc | aucers) |
| Kons Comonnity | As per uneclive 2011/03/EU | | |

*) parameterizable
1) Measuring ranges depend on Memosens sensor
2) As per DIN EN 60746-1, with rated operating conditions
3) ± 1 digit
4) Plus sensor errors
5) c = Cell constant

HandyLab Pilot - PC-Software

HandyLab Pilot is the ideal interface between Memosens[®] and the PC

- Windows interface for intuitive operation.
- The software is automatically connected to HandyLab 750, 750EX and 780 when the devices are connected to the PC.
- Automatic detection of the measuring devices
- If multiple measuring devices are connected, the user can select which one they would like displayed.
- Easy management and evaluation of the measuring data.
- Simple configuration of multiple measuring devices as well as update of the device software if needed.
- Device configurations can be saved and transferred to other devices.
- Set minimum and maximum limit values for the measuring parameters.
- Set difference values.
- Enter your own buffer sets.
- Display device and sensor information (with pH in the shape of a sensor net diagram).
- Save data as CSV files or export into Microsoft[®] Excel.
- Print function.





Order information for HandyLab 7series

| Туре No. | Order No. | Description |
|---------------------|-----------|--|
| HL700AL90pH | 285205110 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with AquaLine 90 pH, Z544, K1A and DIN buffer solutions in vials |
| HL700N1052A | 285205120 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with N1052A and DIN buffer solutions in vials |
| HL700A7781120NMSN | 285205130 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with A7781-120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL700H8281120NMSN | 285205140 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with H8281-120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL 700-PL83120NMSN | 285205150 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with PL83120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL700SL83120NMSN | 285205160 | Set pH-Meter HandyLab 700 for analog and Memosens® pH electrodes with SL83120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL750AL90pH | 285205180 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with AquaLine 90 pH, Z544, K1A and DIN buffer solutions in vials |
| HL750N1052A | 285205190 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with N1052A and DIN buffer solutions in vials |
| HL750A7781120NMSN | 285205200 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with A7781-120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL750H8281120NMSN | 285205210 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with H8281-120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL750PL83120NMSN | 285205220 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with PL83120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL750SL83120NMSN | 285205230 | Set pH-Meter HandyLab 750 for analog and Memosens® pH electrodes with SL83120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL750EXA7781120NMSN | 285205250 | Set pH-Meter HandyLab 750 EX for analog and Memosens® pH electrodes with A7781-120 NMSN, Z544, NMSN1M8EX and DIN buffer solutions in vials |
| HL750EXH8281120NMSN | 285205260 | Set pH-Meter HandyLab 750 EX for analog and Memosens® pH electrodes with J8281-120 NMSN, Z544, NMSN1M8EX and DIN buffer solutions in vials |
| HL750EXPL83120NMSN | 285205270 | Set pH-Meter HandyLab 750 EX for analog and Memosens® pH electrodes with PL83120 NMSN, Z544, NMSN1M8EX and DIN buffer solutions in vials |
| HL750EXSL83120NMSN | 285205280 | Set pH-Meter HandyLab 750 EX for analog and Memosens® pH electrodes with SL83120 NMSN, Z544, NMSN1M8EX and DIN buffer solutions in vials |
| HL780AL90pH | 285205320 | Set pH-Meter HandyLab 780 for analog and Memosens® pH electrodes with AquaLine 90 pH, Z544, K1A and DIN buffer solutions in vials |
| HL780N1052A | 285205330 | Set pH-Meter HandyLab 780 for analog and Memosens® pH electrodes with N1052A and DIN buffer solutions in vials |
| HL780A7781120NMSN | 285205340 | Set pH-Meter HandyLab 780 for analog und Memosens® pH electrodes with A7781-120 NMSN, Z544, NMSN1M8 and DIN-buffer solution in vials |

| Туре No. | Order No. | Description |
|-------------------|-----------|---|
| HL780H8281120NMSN | 285205350 | Set pH-Meter HandyLab 780 for analog and Memosens® pH electrodes with H8281-120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL780PL83120NMSN | 285205360 | Set pH-Meter HandyLab 780 for analog and Memosens® pH electrodes with PL83120 NMSN, Z544, NMSN1M8 and DIN buffer solutions in vials |
| HL780SL83120NMSN | 285205370 | Set pH-Meter HandyLab 780 for analog and Memosens® pH electrodes with SL83120 NMSN, Z544, NMSN1M8 and DIN buffer solution in vials |

Order information for electrodes

| Туре No. | Order No. | Description |
|----------------|-----------|--|
| A7781-120 NMSN | 285114765 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, 3 x ceramic diaphragm, gel electrolyte, Silamid® reference system, ball membrane, A glass, Memosens® threaded plug head, length 120 mm, 12 mm Ø, -5+80 °C, 014 pH |
| AquaLine 90 pH | 285113176 | Plastic shaft, fiber diaphragm, gel electrolyte, Silamid® reference system, cylinder membrane, A glass, threaded plug head, length 120 mm, 12 mm Ø, -5+80 °C, 213 pH |
| H8281-120 NMSN | 285111020 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, KPG ring gap diaphragm, Referid® electrolyte, Silamid® reference system, spherical membrane, A glass, Memosens® threaded plug head, length 120 mm, 12 mm Ø, -5+80 °C, 213 pH |
| N 1052 A | 1054512 | Glass shaft, platinum diaphragm, KCI 3 mol/l electrolyte, Silamid® reference system, temperature sensor Pt 1000, ball membrane, A glass, DIN- +4 mm banana plug, length 120 mm, 12 mm Ø, -5+100 °C, 000.14 pH |
| PL 83-120 NMSN | 285113495 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, Silamid® reference system, 2-hole diaphragm, DuraLid reference system, ball membrane, H glass, Memosens® threaded plug head, length 120 mm, 12 mm Ø, 014, 0130 °C |
| SL 83-120 NMSN | 285114025 | Low-maintenance pH combination electrode with NTC 30K, hot steam sterilizable, CIP and SIP capable, glass shaft, Silamid® reference system, ceramic diaphragm, pressure applied RheoLid reference system, ball membrane, S glass, Memosens® threaded plug head, length 120 mm, 12 mm Ø, 014, 0140 °C |

Order information for accessories

| Туре No. | Order No. | Description |
|-----------|-----------|--|
| NMSN1M8 | 285205380 | Plug/cable combination 1.5 m for digital sensors with Memosens threaded plug head for HandyLab 7series |
| NMSN1M8EX | 285205290 | Plug/cable combination ATEX 1.5 m for digital sensors with Memosens threaded plug head for HandyLab 7series |
| NMSN3M8EX | 285205300 | Plug/cable combination ATEX 3.28 yd for digital sensors with Memosens threaded plug head for HandyLab 7series |
| Z540 | 285205470 | Li-Ion batteries (can only be charged via USB for HandyLab 750 and 780) |
| Z541 | 285205480 | Sensor case (5 ea.) for the liquid-tight storage for HandyLab 7series |
| Z542 | 285205490 | Durable field case for storage and transport for HandyLab 7series and sensor |
| Z543 | 285205500 | Temperature sensor Pt1000 for HandyLab 7series |
| Z544 | 285205510 | 3 adapters for storage of sensors with 12 mm diameter and Pg13.5 thread in the HandyLab 7series case Z541 |

Memosens[®] - Sensors

Mobile multi-parameter Memosens® world with data storage

Our sensors, e.g. the H8281HD, have been proven for decades and make up the basis of our Memosens[®] sensors. Analog measurements are processed in the Memosens[®] plug head and converted to interference-resistant digital signals. When the sensors are connected to the measuring device, they automatically identify themselves with their serial number and part number and process the calibration data stored in the sensor.

Advantages producing reliable measurements

- When connected with HandyLab 7series, there are several advantages for the applications of e.g. samplers in the process. These often contain a multitude of sensors to cover your measurement tasks. Compared to conventional systems, there is no need to calibrate the sensor for every replacement as the calibration data is ow saved in the sensor itself.
- The possibility to distinguish between sensors of the same type by their serial numbers allows for easy allocation and documentation of electronically recorded and saved measurement results.

- In addition to the facilitation of the documentation tasks, the uncertainties during the calibration, such as "hidden serial numbers" etc. are drastically reduced pesky and time-consuming tasks which hinder the operation and favor errors.
- The sensors used to measure pH and ORP (Redox) deliver very small electrical voltages, which can only be measured safely over a limited length (a few meters) of cable when using analog models. Especially outside the laboratory, the measuring points are often difficult to access. Here, the Memosens[®] concept assists by having the calibration data stored in the sensor and via a digital signal transmission.



Memosens® - Sensors - Order Information

| Туре No. | Order No. | Description |
|------------------------|-----------|--|
| A7781-120 NMSN | 285114765 | Low maintenance pH combination electrode with NTC 30K, glass shaft, 3 x ceramic junction, gel electrolyte, Silamid®-reference system, sphere membrane, A-glass, Memosens® screw plug head, length 120 mm, 12 mm Ø, -5+80 °C, 014 pH ¹ |
| A7781-225 NMSN | 285114770 | Low maintenance pH combination electrode with NTC 30K, glass shaft, 3 x ceramic junction, gel electrolyte, Silamid®-reference system, sphere membrane, A-glass, Memosens® screw plug head, length 225 mm, 12 mm Ø, -5+80 °C, 014 pH ¹ |
| FL A 93-120 MF NMSN | 285118180 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 200, sphere membrane, A-glass, Memosens® screw plug head, length 120 mm, 12 mm \emptyset , -30+100°C, 014 pH |
| FL S 93-120 MF NMSN | 285118200 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 300, sphere membrane, S-glass, Memosens® screw plug head, length 120 mm, 12 mm \emptyset , +10+135°C, 014 pH |
| FL A 93-225 MF NMSN | 285118185 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 200, sphere membrane, A-glass, Memosens® screw plug head, length 225 mm, 12 mm \emptyset , -30+100°C, 014 pH |
| FL A 93-280 MF NMSN | 285118190 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 200, sphere membrane, A-glass, Memosens® screw plug head, length 280 mm, 12 mm \emptyset , -30+100°C, 014 pH |
| FL A 93-380 MF NMSN | 285118195 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 200, sphere membrane, A-glass, Memosens® screw plug head, length 380 mm, 12 mm \emptyset , -30+100°C, 014 pH |
| FL S 93-225 MF NMSN | 285118210 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 300, sphere membrane, S-glass, Memosens® screw plug head, length 225 mm, 12 mm \emptyset , +10+135°C, 014 pH |
| FL S 93-280 MF NMSN | 285118220 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 300, sphere membrane, S-glass, Memosens® screw plug head, length 280 mm, 12 mm \emptyset , +10+135°C, 014 pH |
| FL S 93-380 MF NMSN | 285118230 | Liquid electrolyt comb. electrode with NTC 30K and integrated plastic adapter for hose connection incl. Pg13.5 thread, platinum junction, electrolyte L 300, sphere membrane, S-glass, Memosens® screw plug head, length 380 mm, 12 mm \emptyset , +10+135°C, 014 pH |
| H 8281-120 NMSN | 285111020 | Low maintenance pH combination electrode, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, dome membrane, H-glass, Memosens® screw plug head, length 120 mm, 12 mm Ø, 0+100 °C, 213 pH ¹ |
| H 8281-225 NMSN | 285111040 | Low maintenance pH combination electrode, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, dome membrane, H-glass, Memosens® screw plug head, length 225 mm, 12 mm Ø, 0+100 °C, 213 pH ¹ |
| H 8281-325 NMSN | 285111050 | Low maintenance pH combination electrode, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, dome membrane, H-glass, Memosens® screw plug head, length 325 mm, 12 mm Ø, 0+100 °C, 213 pH ¹ |
| H 8281-360 NMSN | 285111030 | Low maintenance pH combination electrode, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, dome membrane, H-glass, Memosens® screw plug head, length 360 mm, 12 mm Ø, 0+100 °C, 213 pH ¹ |
| H 8281-425 NMSN | 285111060 | Low maintenance pH combination electrode, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, dome membrane, H-glass, Memosens® screw plug head, length 425 mm, 12 mm Ø, 0+100 °C, 213 pH ¹ |
| PL 83-120 NMSN | 285113495 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sphere membrane, H glass, Memosens® screw plug head, length 120 mm, 12 mm Ø, 014 pH, 0130 °C 1 |

| Туре No. | Order No. | Description |
|---------------------|-----------|---|
| PL 83-225 NMSN | 285113505 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sphere membrane, H glass, Memosens® screw plug head, length 225 mm, 12 mm Ø, 014 pH, 0130 °C 1 |
| PL 83-325 NMSN | 285113515 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sphere membrane, H glass, Memosens® screw plug head, length 325 mm, 12 mm Ø, 014 pH, 0130 °C 1 |
| PL 83-360 NMSN | 285113525 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sphere membrane, H glass, Memosens® screw plug head, length 360 mm, 12 mm Ø, 014 pH, 0130 °C 1 |
| PL 83-425 NMSN | 285113535 | Low-maintenance pH combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sphere membrane, H glass, Memosens® screw plug head, length 425 mm, 12 mm Ø, 014 pH, 0130 °C 1 |
| PL 89-120 NMSN | 285113565 | Low-maintenance metal combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sensor platinum disk, Memosens® screw plug head, length 120 mm, 12 mm Ø, 014 pH, 0130 °C ¹ |
| PL 89-225 NMSN | 285113575 | Low-maintenance metal combination electrode with NTC 30K, glass shaft, SILAMID® reference, 2 hole junction, DuraLid reference system, sensor platinum disk, Memosens® screw plug head, length 225 mm, 12 mm Ø, 014 pH, 0130 °C 1 |
| Pt 8281-120 NMSN | 285111070 | Low maintenance redox-combination electrode with NTC 30K, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, sensor platinum disk 6 mm Ø, Memosens® screw plug head, length 120 mm, 12 mm Ø, -5+100 °C ¹ |
| Pt 8281-225 NMSN | 285111110 | Low maintenance redox-combination electrode with NTC 30K, glass shaft, KPG annular gap junction, Referid®-electrolyte, Silamid®-reference system, sensor platinum disk 6 mm Ø, Memosens® screw plug head, length 225 mm, 12 mm Ø, $-5+100$ °C ¹ |
| S26250 NMSN | 285128350 | Liquid electrolyte electrode with NTC 30K and KCl vessel 130 mm, Glass shaft, ceramic junction, electrolyte L 310, Ag/AgCl-reference system, cylinder membrane, S glass, length 250 mm, 12 mm Ø, 014 pH, +10+135 °C ¹ |
| SL 83-120 NMSN | 285114025 | Low-maint. pH comb. electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic diaphr., press. RheoLid ref. system, sphere membrane, S glass, Memosens® screw plug head, length 120 mm, 12 mm \emptyset , 014 pH, 0140 °C ¹ |
| SL 83-225 NMSN | 285114035 | Low-maint. pH comb. electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic diaphr., press. RheoLid ref. system, sphere membrane, S glass, Memosens® screw plug head, length 225 mm, 12 mm \emptyset , 014 pH, 0140 °C ¹ |
| SL 83-325 NMSN | 285114045 | Low-maint. pH comb. electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic diaphr., press. RheoLid ref. system, sphere membrane, S glass, Memosens® screw plug head, length 325 mm, 12 mm \emptyset , 014 pH, 0140 °C ¹ |
| SL 83-360 NMSN | 285114055 | Low-maint. pH comb. electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic diaphr., press. RheoLid ref. system, sphere membrane, S glass, Memosens® screw plug head, length 360 mm, 12 mm \emptyset , 014 pH, 0140 °C ¹ |
| SL 83-425 NMSN | 285114065 | Low-maint. pH comb. electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic diaphr., press. RheoLid ref. system, sphere membrane, S glass, Memosens® screw plug head, length 425 mm, 12 mm \emptyset , 014 pH, 0140 °C ¹ |
| SL 89-120 NMSN | 285114075 | Low-maintenance metal comb. Electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic junction, press. RheoLid ref. system, sensor platinum disk, Memosens® screw plug head, length 120 mm, 12 mm Ø, 0140 °C |
| SL 89-225 NMSN | 285114085 | Low-maintenance metal comb. Electrode with NTC 30K, hot steam sterilizable, CIP and SIP able, glass shaft, SILAMID® ref., ceramic junction, press. RheoLid ref. system, sensor platinum disk, Memosens® screw plug head, length 225 mm, 12 mm Ø, 0140 °C |

All electrodes have in common:

Temperature sensor NTC30kΩ
Glass shaft, diameter 12mm

What can Xylem do for you?

We're 12,700 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to www.xyleminc.com



a xylem brand

SI Analytics GmbH

Hattenbergstr. 10 55122 Mainz Germany

Phone:+49 6131 66 5111Fax:+49 6131 66 5001E-Mail:si-analytics@xyleminc.comInternet:www.si-analytics.com

presented by

SI Analytics is a trademark of Xylem Inc. or one of its subsidiaries. © 2014 Xylem, Inc. **980 080US** Version 06/2014