

Tergeo plasma cleaners

Versatile tabletop plasma cleaner for R&D and low-volume production



Customers around the world



HARVARD
UNIVERSITY



Massachusetts
Institute of
Technology



Berkeley
UNIVERSITY OF CALIFORNIA



THE
UNIVERSITY OF
BRITISH
COLUMBIA
Canada



NUS
National University
of Singapore
Singapore



Australian
National
University

Australia



Germany

MAX PLANCK INSTITUT
FÜR MIKROSTRUKTURPHYSIK

United Kingdom
UNIVERSITY OF LEEDS



France



清华大学
Tsinghua University
China



Weizmann Institute of Science
Israel

جامعة الملك عبد الله
للعلوم والتقنية
King Abdullah University of
Science and Technology

Saudi Arabia



UNIVERSIDADE FEDERAL
DO RIO GRANDE DO SUL
Brazil



고려대학교
KOREA UNIVERSITY

South Korea



TOKYO METROPOLITAN UNIVERSITY
首都大学東京

Japan

Certified by TÜV Rheinland and UL Japan

IEC 61010-1:2010 certificate for **European countries**

UL CAN/CSA 61010-1:2012 certificate for **North America**

Rf emission conformity certificate from **UL Japan Inc.**


		Ref. Certif. No. US-TUVR-012083
IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB SCHEME		
CB TEST CERTIFICATE		
Product	Plasma Cleaner	
Name and address of the applicant	PIE Scientific LLC 3209 Whipple Rd Union City CA 94587, USA	
Name and address of the manufacturer	PIE Scientific LLC 3209 Whipple Rd Union City CA 94587, USA	
Name and address of the factory	PIE Scientific LLC 3209 Whipple Rd Union City CA 94587, USA	
Ratings and principal characteristics	AC 100-230V; 50/60Hz; 300W; Class I	
Trademark (if any)	PIE Scientific	
Customer's Testing Facility (CTF) Stage used		
Model / Type Ref.	Tergeo, Tergeo Plus, Tergeo-EM	
Additional information (if necessary may also be reported on page 2)	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> CE marked </div>	
A sample of the product was tested and found to be in conformity with	IEC 61010-1:2010 for national deviations see test report	
As shown in the Test Report Ref. No. which forms part of this Certificate	31973182 001	
This CB Test Certificate is issued by the National Certification Body		
	TÜV Rheinland of North America, Inc. 12 Commerce Road Newtown, CT 06470 - USA	Date: 26.03.2020 Signature:

<h2>Certificate</h2>	
Certificate no.	CU 72200728 01
License Holder:	PIE Scientific LLC 3209 Whipple Rd Union City CA 94587 USA
Manufacturing Plant:	PIE Scientific LLC 3209 Whipple Rd Union City CA 94587 USA
Test report no.:	USA- 31973182 003
Client Reference:	
Tested to:	UL 61010-1:2012 R4.16 CAN/CSA-C22.2 NO. 61010-1-12 + GI1 + GI2 (R2017)
Certified Product:	Plasma Cleaner
License Fee - Units	
Model Designation:	Tergeo, Tergeo Plus, Tergeo-EM 7
Rated Voltage:	AC 100-230V, 50/60Hz
Rated Power:	300W
Protection Class:	I
<div style="border: 1px solid black; padding: 10px; display: inline-block;"> TUV Rheinland marked for OSHA NRTL compliant in North America </div>	
Licensed Test mark:	
Date of Issue (day/mo/yr)	03/08/2020
TÜV Rheinland of North America, Inc., 12 Commerce Road, Newtown, CT 06470, Tel: (203) 426-0888 Fax: (203) 426-4000	


	Test report No. : 12402432S-B Page : 1 of 18 Issued date : August 28, 2018
RADIO TEST REPORT	
Test Report No.: 12402432S-B	
Applicant	: PIE Scientific LLC
Type of Equipment	: Plasma cleaner
Model No.	: Tergeo-EM
Test standard	: Rules for Regulating Radio Equipment, Article 65-3 (1) : Public Notice of Ministry of Internal Affairs and Communications, No. 207
Test result	: Complied
1. This test report shall not be reproduced in full or partial, without the written approval of UL Japan, Inc. 2. The results in this report apply only to the sample tested. 3. This sample tested is in compliance with the limits of the above standard. 4. The test results in this test report are traceable to the national or international standards. 5. This test report covers Radio technical requirements. It does not cover administrative issues such as Manual or non-Radio test related Requirements. (if applicable) 6. The all test items in this test report are conducted by UL Japan, Inc. Shonan EMC Lab.	
Date of test:	December 4, 2017 and August 10, 2018
Representative test engineer	 Akihiro Oda Engineer Consumer Technology Division
Approved by:	 Toyokazu Inamura Leader Consumer Technology Division
UL Japan, Inc. Shonan EMC Lab. 1-22-3 Megurogaoka, Hiratsuka-shi, Kanagawa-ken, 259-1220 JAPAN Telephone : +81 463 50 6400 Facsimile : +81 463 50 6401	

Global sales coverage

Made in San Francisco
Bay Area in U.S.A.



Sold to over 30 countries and regions



North America: The United States, and Canada

Europe: United Kingdom, Germany, France, Italy, Netherlands, Austria, Sweden, Switzerland, Luxembourg, Czech Republic, Poland, Hungary, Croatia, and Monaco

Asia: Japan, S. Korea, China, Hong Kong, Taiwan, Singapore, Malaysia, India, and Thailand

Middle East: Israel and Saudi Arabia


Oceania: Australia, and New Zealand

South America: Brazil, and Peru

Key Spec

Items	Feature and Spec
RF power supply	75-watt or 150-watt for Tergeo; 150-Watt or 500-Watt for Tergeo-plus; can be adjusted at 1-watt interval; 13.56MHz high frequency RF power supply; automatic impedance matching; Continuous and pulsed plasma (pulse ratio 100% to <1%) .
Sample chamber	Cylindrical quartz chamber (ID:110mm, OD:120mm, Depth 280mm for basic Tergeo; ID160mm, OD170mm, Depth 280mm for Tergeo plus). One rectangular quartz shelf included.
Plasma discharge mode	Capacitively coupled discharge for the direct plasma source. Inductively coupled plasma source as the remote plasma source.
Electrode placement	External rf electrode design to reduce ion sputtering on the metal electrodes.
Source design	Dual plasma sources in one system . Both the traditional direct mode plasma processing and the gentle downstream mode plasma processing in one system
Plasma diagnostics	Realtime plasma intensity sensor for quantitative plasma status monitoring
Gas delivery	One solenoid valve controlled gas input port for chamber venting/purging purpose; two or three additional MFC regulated gas input ports for process gas delivery. All the gas input ports use the standard 1/4" Swagelok tube compression fitting.
Pressure sensor	Premium full range corrosion resistant (platinum filament) pressure sensor.
Pump	Basic oxygen service oil pump (ultimate pressure <10mTorr) or oil-free dry pump from Edwards Vacuum LLC. Pump can be controlled by the plasma system.
System control	7 inch LCD touchscreen. 20-recipe library. Fully automatic operation. 3-recipe job sequence mode.
Enclosure dimension	W450mm X H250mm X D430mm for basic Tergeo, W500mm X H300mm X D430mm for Tergeo-plus
AC input	Universal 100V~230V AC input 50/60Hz

Better uniformity and lower contamination



Inductive coil antenna design in some low-cost plasma cleaner results in bad uniformity

Better uniformity. Plasma discharge technology in Tergeo plasma cleaner originates from the research carried out in the Plasma & Ion Source Technology at the Lawrence Berkeley National Laboratory. The difference in plasma color is because of the different process gas. Tergeo plasma cleaner clearly has much better uniformity. Bad uniformity means bad process control. Tergeo system also have larger usable volume for the same chamber diameter.


Lower contamination: Electrodes are placed outside of the quartz tube in Tergeo plasma cleaner. Ions can't reach the metal electrodes. If high voltage metal electrode is placed inside the plasma chamber, ions will be accelerated to high energy and sputter metal out of the electrodes. Metal will then deposit onto the samples and cause metal contamination issue. Metal sputtering contamination issue is especially severe if KHz rf power supply is used to generate the plasma because ions can be accelerated by KHz electric field.

Tergeo vs other plasma cleaners

High voltage rf electrodes are placed outside the plasma chamber in Tergeo system. It results in better plasma uniformity, less contamination and more usable space


Tergeo

Direct mode




Better uniformity, no metal contamination

Downstream mode





Gentle ion-free chemical reaction

Competitors



High voltage rf electrode in the plasma attracts high energy ions and cause metal sputtering contamination on the samples

Direct vs downstream cleaning modes



Direct/immersion mode plasma cleaning


Plasma is generated in sample chamber
Samples are immersed in plasma
Samples are subject to chemical reaction with radicals and energetic ion sputtering

Remote/downstream mode plasma cleaning

Plasma is generated outside sample chamber.
Samples are not immersed in energetic plasma
Only gentle chemical reaction takes place on sample surface.
No energetic ion sputtering. No sample damage due to excessive heating, charging and ion bombardment.


Proprietary plasma sensor technology, quantitative data for repeatable process

Direct/immersion cleaning mode



Plasma generated in the sample chamber. Plasma intensity reading is 206


Remote/downstream cleaning mode



Plasma generated in the remote plasma source. Plasma intensity reading is 170

Take the guesswork out of plasma cleaning. Plasma sensor tells the user exactly how strong the plasma is. User can use the quantitative plasma intensity measurement as feedback to adjust gas flow rate and rf power for desired cleaning speed. User doesn't need to be plasma expert to set up right cleaning recipe for different samples.

Continuous and pulsed plasma



Tergeo plasma system can generate **continuous** or **pulsed** plasma.

Tergeo plasma system not only adjusts rf power but also adjusts the pulse duty ratio. Pulse mode further extends the dynamic range of the plasma intensity. The pulse duty ratio can be adjusted from 100% (continuous) to less than 1% (two order of magnitude). The plasma intensity varies more linearly with the pulse duty ratio compared with adjusting the rf power wattage. Pulse mode can create an extremely weak plasma that is almost invisible to human eyes. Our proprietary plasma sensor can still measure the plasma intensity quantitatively even for the plasma that is almost invisible to human eye.

Pulsed mode can generate weak yet repeatable plasma for fragile samples like 2D materials, TEM holey carbon grids, and other samples that can be easily damaged by strong plasma.

KHz, MHz, and GHz plasma

Why MHz is the best frequency to generate the plasma?

- The frequency of kHz is so slow that both electrons and ions in the plasma can be accelerated by the kHz electric field. Therefore, ions and electrons can be accelerated to very high energy in each cycle. They will eventually hit the chamber wall and electrodes and get lost from the plasma. Therefore, KHz plasma is usually very weak. **KHz rf power supply and plasma systems are generally much cheaper than the 13.56MHz rf plasma system due to the lower efficiency.**
- The frequency of the MHz electromagnetic wave is slower enough to accelerate the electrons. But ions are too heavy and too slow to respond to the MHz electric field. Ions are thermalized in the MHz plasma system. So they can survive in the plasma much longer. MHz is the sweet spot of the frequency to generate high-intensity uniform plasma.
- When the frequency increases to GHz (microwave), energy can only be transferred to the plasma through a plasma wave. Microwave tends to create multiple nodes in the metal chamber, which generates localized hot spots and results in bad uniformity. This is why dishes need to be rotated in the microwave oven to improve uniformity. Due to the bad uniformity, microwave plasma usually works in downstream mode.

Intuitive design and easy to use

- Two different user interface: LCD touchscreen user interface and PC remote control user interface.
- Intuitive user interface. Easy to use. No extensive training required.
- Fully automatic operation. Supports 20 recipes. One button to start. Repeatable and reliable results.
- Recipe lock/unlock function to prevent unintended modification of the optimized recipes. Essential features for production runs and for shared facility.
- Smart system self-diagnostic. Notify user the operation tips, system status and problems on the LCD screen or PC user interface.

Supports
total 20
recipes

#	Recipe Name	Mode Selection	Power (Watt)	Duty Rate	Gas 1 (sccm)	Gas 2 (sccm)	Gas 3 (sccm)	Duration (M:S)	Purging Gas	Start (1-9)
1	Pure Argon	Immersion	150	255	3.0	0.0	0.0	2 : 0	No	Yes
2	Pure Oxygen	Immersion	150	255	0.0	2.0	0.0	2 : 0	No	Yes
3	Low power mix	Immersion	25	255	5.0	2.0	0.0	2 : 0	No	Yes
4	Medium pow. mix	Immersion	75	255	5.0	2.0	0.0	2 : 0	No	Yes
5	High pow. mix	Immersion	150	255	2.0	2.0	0.0	2 : 0	No	Yes
6	Remote mix MP	Remote	75	255	5.0	2.0	0.0	2 : 0	No	No
7	Remote MIX HP	Remote	150	255	20.0	3.0	0.0	2 : 0	No	No
8	Remote O2 HP	Remote	150	255	0.0	2.0	0.0	2 : 0	No	No
9	Remote Ar LP	Remote	25	255	10.0	0.0	0.0	2 : 0	No	No
10	Low power O2	Immersion	25	255	0.0	1.0	0.0	2 : 0	No	Yes


Recipe can
be locked for
production

The screenshot shows a recipe configuration window with the following fields and values:

- Number: 3
- Name: O2 PDMS/Glass 2
- Mode: Direct
- RF Power (W): 40
- Pulse ratio (N/255): 255
- Gas 1 (sccm): 0.0
- Gas 2 (sccm): 40.0
- Gas 3 (sccm): 0.0
- Duration(M:S): 0 : 20
- Flush valve: No purge

A red warning box at the bottom states: "This recipe has been locked by the administrator. Click the left button to unlock it if you are the administrator". A "Locked" indicator is visible at the bottom left, and a "Cancel" button is at the bottom right.

Powder sample support




Slow pumping and slow venting features to prevent the powder and light samples from being blown away by strong air flow during pumping and venting. Slow pumping and slow venting feature can be enabled in the touchscreen controller



Spare stepper motor driver to control an external stirring motor for powder samples, improve process uniformity for powder samples.

System log and run record

1. The capability to keep a run log is essential for some bio-medical, pharmaceutical and defense industries.
2. The PC control software provided by PIE Scientific LLC can save the system parameters and run logs in to a CSV file. The CSV file can be processed in excel.
3. PC communicates with Tergeo through an RS232 serial port. A regular USB port on PC can be easily converted to an RS232 port using an USB-RS232 converter. Then the user can simply connect the PC and Tergeo with a D-sub 9 cable.




Data log setting


Process monitoring for production

- Monitor the process through the PC control software.
- Make sure the process is stable for day to day production.
- Alert the operator if the process drifts out of the upper and lower control limit.
- Alert can be set on rf power, chamber pressure, gas flow rate on each channel, plasma intensity reading.
- Process monitor limit can be set for each recipe.

Sample warning message on chamber pressure



Sample warning message on plasma intensity




How do we support customers all over the world?

More than 70% of our sales are international export. Our customers cover North America (U.S.A and Canada), U.K., E.U (such as Germany, France, Sweden, Netherlands, etc.), Asia (Japan, S. Korea, China, Hong Kong, Taiwan, Singapore, Malaysia, Thailand, and India), Australia, New Zealand, and Israel. Reliability is the No.1 priority in our system design. Additionally, our system is designed to be easy to support globally. Here are several advanced features designed for easy to support and service.


- 1) Smart system self-diagnostics capability. The system constantly monitors the status of each sub-components. If it detects an issue, it will display a message on the LCD touchscreen display and prompt user what to do.
- 2) Remote control of the Tergeo plasma system from our U.S. headquarter. Tergeo system can be operated from an external computer. We can log onto the computer at the customer site from the United States to operate the machine and figure out any issues. Skype is one of the remote access platforms.
- 3) PC control software for Tergeo keeps a system status log. The customer can just send us the log file for us to figure out any system issue.
- 4) Remote video support through Skype or Facetime. We can view the LCD touchscreen of the Tergeo system and figure out any issue. We can also guide the user to service and repair the system through remote video support.



Smart system self-diagnostics



Fail to pump down. Wipe off the dust on the door o-ring! Perform helium leak check!



Tip: Powder samples? Enable slow pumping and slow venting in the system screen.

System self-diagnostics:

The system constantly monitors the status of each sub-component. Once it detects an issue, it will display a message on the LCD touchscreen display and prompt user what to do.

For example: one of the previously most reported issue from our customer is the dust particles on the door o-ring that cause vacuum leak. It can be easily resolved by wiping off the dust with a lint free wipe. The MKS pressure sensor used in Tergeo is sensitive to helium gas. Spraying helium gas around the vacuum joint can detect the leak point.

System operation guide:

When the system is in the idle state, the status display window will cycle through pre-stored tip messages and help the user to get familiar with the system faster.

System can also detect how the system is being operated (such as manual mode or recipe mode). It will provide suggestion on how to operate the system in the optimal method based on how the system is being used by the operator.

Remote support methods

Method 1: remotely access the machine from the U.S. headquarter in the Silicon Valley



Customer site: Control the Tergeo with a PC



U.S. headquarter: remotely control the Tergeo through skype or Teamviewer remote screen sharing.

Method 2: remote video support through Skype or Facetime



View the machine in real-time to guide the user to troubleshoot the issue.

Method 3: Upload the data log file to the U.S. headquarter for analysis

Log File Directory	C:\temp\2018_09_10_180920TG97.csv	Change Directory
Data Logging	Data logging Active state: <input type="checkbox"/>	1000
<input checked="" type="checkbox"/> Enabled	interval (ms) Idle state: <input type="checkbox"/>	30000
Need support on system issue, run a recipe, send the recent log file to support@piescientific.com with a detailed description of the issue.		



Application guide

Photoresist ashing and descum

Plasma emission intensity reading



Oxygen plasma photoresist ashing speed.


Tergeo-plus, 150W rf power, 4-inch wafer, Shipley 1813 positive resist, cured

O2 flow rate (sccm)	3	5	5	10	20	30
Argon flow rate (sccm)	0	0	5	0	0	0
Chamber pressure (mTorr)	65	92	123	157	250	336
Plasma intensity reading	298	263	820	210	157	144
Etching speed (nm/min)	186.57	231.17	183.00	160.20	97.20	34.60

[Gentle downstream/remote mode is available as an upgrade option](#)

Tergeo plasma system can measure plasma emission light intensity. The plasma emission light can vary a lot with different gas species. During O₂ plasma photoresist ashing process, it can generate CO, CO₂, H₂O type byproduct, which can change the plasma emission intensity a lot. Therefore, it's possible to monitor the progress of the ashing process by monitoring the plasma emission intensity.



Shipley 1813 photoresist ashing curve in 150W Tergeo-plus



Wire bonding and flip-chip bonding in IC/MEMS packaging

Remove organic, oxide and fluoride surface contamination, render surface hydrophilic, improve bonding strength and packaging reliability

Silver wire bond pull strength test before and after argon plasma cleaning



Tarnished silver pad looks dull before plasma cleaning. The pad looks significantly brighter and shinier after 3-minute argon plasma cleaning.

YouTube video compares chip packaging before and after argon plasma cleaning:
<https://www.youtube.com/watch?v=axQGydpqUAK>

2-D materials research

Applications:

- Substrate cleaning before growing 2D materials.
- Plasma functionalization without damage the 2D materials using our unique downstream and pulsed plasma.
- 2D materials etching and removal.

Unique features on Tergeo for 2D materials research

1. Fast direct plasma etching mode and gentle downstream plasma cleaning mode in one system.
2. Continuous and pulsed plasma for high intensity and extremely weak plasma, wide dynamic range on the plasma intensity.
3. Direct/downstream and continuous/pulsed mode together can change the etching speed by more than 3 orders of magnitude.
4. Plasma emission intensity sensor for quantitative measurement and control of the plasma intensity.

PDMS and Glass bonding

Bond PDMS/Glass immediately, No need to cure the device in an oven





Figure 1



A layer of PDMS remained on the glass after the top part of the PDMS broke off.


Figure 2

Figure 1, After plasma processing, PDMS bonds to glass immediately, trapped gas at the interface is automatically squeezed out. **No need to cure the PDMS/Glass bond in an oven for a couple of hours.**

Figure 2, One minute after bonding, the user tries to break the PDMS/Glass bond. The bond is so strong that the PDMS/Glass bonding interface remains intact even after the top part of the PDMS was pulled off.

Surface activation for bio-polymer/plastic before coating


Typical process gas: room air (N₂+O₂), pure O₂, argon, or argon+O₂
Processing time: 20 seconds to 5 minutes




Before cleaning,
Water contact angle $>80^\circ$



Right after 50W, 60s
plasma cleaning.
Water contact angle is
too low to be measured.




3 hours after plasma
cleaning.
Water contact angle $=12^\circ$



63 hours after plasma
cleaning.
Water contact angle $=21^\circ$

Glass, wafer, and well plate cleaning

Glass slides




Glass slide surface before cleaning



Glass slide surface after cleaning

Wafer



Wafer surface before cleaning




Wafer surface after cleaning

Well plate



Well plate in ambient air plasma on a quartz holder



Upside down well plate in oxygen plasma


Activate surface and remove organic contamination.

Rendering the glass slides to super hydrophilic, then to strongly hydrophobic

Hydrophilic surface treatment



Original glass slides




Right after 1-minute air plasma treatment

Hydrophobic surface treatment





Right after proprietary hydrophobic plasma surface treatment



15 hours after hydrophobic surface treatment

Medical implant cleaning and activation

Dental implants before plasma processing



Video showing the water contact angle before and after the pure O₂ plasma processing:

Before plasma processing:


https://www.piescientific.com/Customer_demo_results/Medical_Implant/75W_pure_oxygen/Before.MOV

After argon plasma processing:


https://www.piescientific.com/Customer_demo_results/Medical_Implant/75W_pure_oxygen/After.MOV

SEM/TEM sample cleaning


Before cleaning




After 1-minute cleaning




Hydrophobic carbon grid



Hydrophilic carbon grid




Ultra-thin TEM carbon grid
After plasma activation




Decapsulation for SEM failure analysis

1. Remove parylene conformal coating in PCB failure analysis with oxygen plasma
2. Remove polyimide capping layer with oxygen
3. Remove nitride and oxygen capping layer with CF_4+O_2 plasma

Remove 5um parylene coating on a glass slide sample in 20 minutes



Kapton tape (1 mil (25um) polyimide film over 3 mil silicone adhesive)
25um polyimide film is mostly remove after 3-4 hours of oxygen plasma. Only the silicone adhesive layer remained.




Plasma ashing to remove organic materials



Paper in Ar+O₂ plasma



After 10 hours of plasma ashing



Only inorganic matters are left for analysis

Remove organic materials to reveal the inorganic matters. For example, analyzing asbestos in paint, content of the inorganic matters in coal powder.

- The change of the plasma color is due to the change of the ashing by-products CO_x and H₂O. Our unique plasma intensity sensor reading can pick up the change of the plasma emission.
- Slow venting option to prevent the ash powder from being blown away during venting.