




**EDEN LAB**  
LABORATORY SUPPLIES



Authorised Distributor  
For



[contact@edenlab.com.au](mailto:contact@edenlab.com.au)  
[www.edenlab.com.au](http://www.edenlab.com.au)



MDS1  
Middle School Microscope



## **Magnification Options**

Supplied with a Widefield 10x eyepiece, along with 4x, 10x, 40xr (retractable) objectives to provide a range of magnification options: 40x, 100x, 400x.

## **Ease of Use**

Coarse and Fine focus knobs make it easy to obtain a crisp and clearly focused image. Rack stop protects objectives and specimen slides from damage. Head rotates 360°.

Specimen circular stage slides forward, backward and left-right. Locked-on spring-activated stage clips hold samples in place.

Rugged, all-metal frame has enamel coating for protection. Carrying handle protects microscope from damage during transport.

## **Cordless or Cordless LED Illumination**

Cool, bright LED Illumination. Includes rechargeable batteries and recharger. Use corded or cordless. CE certified power cord. 100~240V. Rheostat control for illumination.

## **You're Covered!**

Includes dust cover and 5-year warranty.



MDS1 Microscope Magnification

Achromat Objective	10x Eyepiece
4x	40x
10x	100x
40x	400x

- Autoclaves
- Balances & Scales
- Bath & Circulators
- Centrifuges
- Chillers/ Freeze Dryers
- Desiccators
- Filtration Systems
- Furnace/ Burners
- General Lab Equipment
- Glovebox
- Heating & Cooling blocks
- Heating mantle/ Temperature controls
- Homogenizers
- Hotplate & Stirrers
- Incubators
- Lab Cabinets & Storage
- Laboratory consumables
- Measurement & Testing  
(pH/ DO/ TDS/ Refractometer/Conductivity)
- Microscope (Research)
  - o Life science Microscope
  - o Geology Microscope
  - o Industrial Inspection Microscope
  - o Material Science Microscope
  - o Microscope Cameras
- Microscope (Educational/ Teaching)
  - o Highschool
  - o College
  - o University
- Ovens
- Pipette
- Plasticwares
- Rangehoods/ Clean Bench, BIO Safety Cabinet
- Recirculating Coolers
- Refrigerator/ Freezers
- Rotary Evaporator
- Rotators
- Shakers/ Mixer/ Roller
- Temperature/ Environment Chambers
- Ultrasonic cleaner
- UV Vis Spectrophotometer
- Vacuum Pump

